




MONOCULAR MICROSCOPE



SIZE AND WEIGHT

	SLV-10M	SLV-10B	SLV-10D
Product Size	18x13x36CM	18x13x36CM	18x13x36CM
Net Weight	2.25KG	3.14KG	3.2KG



Technical Specification	
MODEL No.	SLV-10M
Viewing Head:	Monocular, 45° inclined; 360° rotating.
Eyepiece:	WF10x/18 mm. OPTINAL: Pointer, Micrometre Eyepiece
Nosepiece:	Quadruple ball bearings revolving nosepiece.
Objectives:	Achromatic 4x/0.10, with anti-fungus treatment Achromatic 10x/0.25, with anti-fungus treatment Achromatic 40x/0.65, with anti-fungus treatment Achromatic 100x/0.17, with anti-fungus treatment
Specimen stage:	Double layer Mechanical Stage 110mm X 120mm. / 60mm X 30mm.
Focusing:	Coaxial coarse and fine focusing mechanism with limit stop to prevent the contact between objective and specimen.
Condenser:	Diffusing filter with rotating diaphragm wheel 1.25NA
Illumination:	1W LED, with light intensity control / LED with Battery Back-up.

OPTIONAL: Pointer Attachment, Micrometre/Grid/Grain Attachment, DIGITAL Attachment.



Shop No.1, Basement floor, Best Apartments, Dr. Arunachalam Road, 5th Stage BEML Layout, Rajarajeshwari Nagar, Bengaluru - 560 098