


10.5" Touch Screen Android Pad Digital Camera, 8.0M



10.5" LCD Touch Screen Digital Microscope, 8.0M, Android Pad



- 10.5" Touch Screen Android Pad + 4K 8.0M Digital All-in-One Camera for Microscope
- Turn On & See Micro Image Directly Through Pre-Installed S-Eye Measure Software
- Control by Touch Screen & USB Mouse+Keyboard to Measure No Need Computer
- 4K 8.0M 1/1.8" CMOS Digital Camera Resolution 3840*2158@30FPS Shape Image
- Built-in Android 11 System Support 90 Language, With Memory RAM 2G + ROM 16G





10.5" Android Pad + 8.0M Camera
Sharp Image | All-in-One | Phot & Video





10.5" Touch Screen Android Pad Digital Camera, 8.0M

Spec.	Model	SLV-8MP-Android
	Category	HD All-in-One Android 11 System Tablet Digital Camera
	Sensor	1/1.8" CMOS Sensor
	Effective Pixel	8M, 4K
	Screen	10.5" High-Definition Touch Screen
	Operation System	Android 11
	WIFI	2.4GHz/5GHz Dual-band WIFI, Support WIFI6
	Bluetooth	Bluetooth 5.0
	Output	USB3.0*2, USB2.0*1, HDMI, Gigabit LAN Port
	Pixel size	2.0um*2.0um
	Power	DC-12V/2A
Function	Basic Function	Android Interaction, Installed S-Eye2.0 Measure Software
	Image Adjustment	Brightness, Contrast, Saturation, Color Temperature
	Camera Function	Mirror, Flip, Freeze, Black and White
	UI Interface	Support User-friendly Mouse & Keyboard, And Touch Screen Operation
	Photo	8.0M Image Capture
	Video	1080P Video
	Storage	Still Photo and Video Recording, Picture Preview
	Marking	Point Coordinates, Crosshairs, Coordinate System, Text Annotations
	Length Measurement	Straight Line, Folded Line, Curve Length, Parallel Line Distance, Point Line Distance
	Geometric Measurement	Line Segment Length, Radius Fixed Circle, 2 Point Fixed Circle, 3 Point Fixed Circle, Concentric Circle, Radius Fixed Circle
Package	Geometric Area Measurement	Polygon, Rectangle
	Measurement System	Built-in S-Eye Software
	Box Size	237*169*10.5mm
Package	Weight	0.6kg
	Accessory	DC-12V/2A Power Adapter



📍 Shop No.1, Basement floor, Best Apartments, Dr. Arunachalam Road, 5th Stage BEML Layout,
 Rajarajeshwari Nagar, Bengaluru - 560 098