

SG HORIZONTAL AUTOCLAVE MAKER



Use

The table type sterilizers are safely auto-controlled which are designed for the dental clinics, hospitals, school lab, scientific research institutions and Laboratory and mostly suitable for sterilizing of surgical, dental and eye instruments, glassware, medicine, culture medium and biological dressing, food etc.

Characteristics

- For 4-6 minutes rapidly sterilizing.
- Digital display of working status, touch type key.
- With 3 fixed cycles of sterilizing program and adjustable for user.
- The whole process of water adding, temperaturerising, sterilizing, drying steam discharging controlled automatically.
- Steam-water inner circulation system: no steam discharge, and the environment for sterilizing will be clean and dry.
- Safe protection of water lacking.
- Exhaust the cool air automatically.
- Door safety lock system.
- With four stainless steel sterilizing plates.

The chamber of the sterilizer is made of stainless steel.

Automatically shut off with beep reminding after sterilization.

With drying function

| MODEL | SG-HAC35D | SG-HAC50D |
|----------------------------|--------------------|--------------------|
| Sterilizing chamber volume | 35L Φ300x500mm | 50L Φ340x550mm |
| Working pressure | 0.22Mpa | |
| Working temperature | 134°C | |
| Temperature adjust range | 105-134°C | |
| Timer range | 0-99min | |
| Heat average | ≤±1°C | |
| Power | 2.7KW/AC220V.50Hz | |
| Sterilizing plate | 400x200x30(4piece) | 500x250x30(4piece) |
| Dimension | 730x550x510 | 770x600x540 |
| Exterior package dimension | 810x680x650 | 880x730x700 |
| G.W/N.W | 90/73kg | 105/86kg |

SKU: HAC-XD35D Category: [Table Top Autoclave](#) Tags: [Autoclave](#), [BJ](#), [Table-Top Autoclave](#)



Shop No.1, Basement floor, Best Apartments, Dr. Arunachalam Road, 5th Stage BEML Layout, Rajarajeshwari Nagar, Bengaluru - 560 098