



LABMICRO

8553545575 / 8050259428 

contact@labmicro.in 

www.labmicro.in 



SLV-622 UP/DOWN Holder

SLV-603 Body with
SLV-621-4C Universe Stand



SLV-621-1 UNIVERSAL STAND

SLV-621-1C UNIVERSAL STAND



SLV-621-2 UNIVERSAL STAND

SLV-621-2C UNIVERSAL STAND



SLV-621-3 UNIVERSAL STAND

SLV-621-3C UNIVERSAL STAND



SLV-621-4 UNIVERSAL STAND

SLV-621-4C UNIVERSAL STAND



Shop No.1, Basement floor, Best Apartments, Dr. Arunachalam Road, 5th Stage BEML Layout,
Rajarajeshwari Nagar, Bengaluru - 560 098

	1	2	3	4	1C	2C	3C	4C
Stand	Soleplate 254x254x40mm				U clamp maximum clamped thickness 63.2mm			
Steeve	The measurement to match between the steeve and the aperture is $\Phi 25\text{mm}$							

	1,1C,2,2C	3,3C,4,4C
Post	Height 475mm, $\Phi 37\text{mm}$	Height 550mm, $\Phi 37\text{mm}$
Rail	Flexed length 560-605mm	Double rail length 460mm, $\Phi 19.05\text{mm}$



8 Stereo Microscope Series + 35 Stands

8 Stereo Microscope Series + 35 Stands

The Head Can Be Combined Freely With Different Stands & Holder Arms



SLV-630 35 kinds stand Optional

SLV-630 35 kinds stand optional, Details see choose model according to catalogue no. Like:
SLV-45-B1 Means Binocular head with B1 stand
SLV-45-B1T Means Trinocular head with B1 stand

1



2



A1 Holder

+



B1 Pole Stand

=



B1 Stand

3



A3 Holder

+



STL4 Boom Stand

=



STL4 Stand

SLV-630 Pole Stand Models



SLV-630-B1/B1A1

SLV-630-B1 Stand 200*255*22mm, Pole Height 240mm, Dia.32mm, No Light Source.

SLV-630-B1A1 Stand 200*255*22mm, Pole Height 240mm, Dia.32mm, No Light Source.
With SLV-640-A1 Holder Inner Dia.76mm.



SLV-630-B2/B2A1

SLV-630-B2 Stand 200*255*60mm, Pole Height 240mm, Dia.32mm, Up/Bottom 12V15W Halogen Lamp.

SLV-630-B2A1 Stand 200*255*60mm, Pole Height 240mm, Dia.32mm, Up/Bottom 12V15W Halogen Lamp. With SLV-640-A1 Holder Inner Dia.76mm.



SLV-630-B2LD

SLV-630-B2LDA1

SLV-630-B2LD/B2LDA1

SLV-630-B2LD Stand 200*255*60mm, Pole Height 240mm, Dia.32mm, Up 30 LED & Bottom 92 LED.

SLV-630-B2LDA1 Stand 200*255*60mm, Pole Height 240mm, Dia.32mm, Up 30 LED & Bottom 92 LED. With SLV-640-A1 Holder Inner Dia.76mm.



SLV-630-B6

SLV-630-B6A1

SLV-630-B6/B6A1

SLV-630-B6 Stand 205x275x22mm, Pole Height 240mm, Dia.32mm, No Light Source.

SLV-630-B6A1 Stand 205x275x22mm, Pole Height 240mm, Dia.32mm, No Light Source. With SLV-640-A1 Holder

Inner Dia.76mm.



SLV-630-B8LS

SLV-630-B8LSA1

SLV-630-B8LS/B8LSA1

SLV-630-B8LS Stand 205x275x40mm, Pole Height 240mm, Dia.32mm, Double Pipe LED.

SLV-630-B8LSA1 Stand 205x275x40mm, Pole Height 240mm, Dia.32mm, Double Pipe LED. With SLV-640-A1 Holder Inner Dia.76mm.



SLV-630-B8L

SLV-630-B8LA1

SLV-630-B8L/B8LA1

SLV-630-B8L Stand 205*275*40mm, Pole Height 240mm, Dia.32mm, Up/Bottom 3W LED.

SLV-630-B8LA1 Stand 205*275*40mm, Pole Height 240mm, Dia.32mm, Up/Bottom 3W LED. With SLV-640-A1

Holder Inner Dia.76mm.



SLV-630-J4L/J4LA1

SLV-630-J4L Stand 285*238*40mm, Pole Height 240mm, Dia.32mm, Up/Bottom 3W LED.

SLV-630-J4LA1 Stand 285*238*40mm, Pole Height 240mm, Dia.32mm, Up/Bottom 3W LED. With SLV-640-A1 Holder Inner Dia.76mm.



SLV-630-J4LD/J4LDA1

SLV-630-J4LD Stand 285*238*40mm, Pole Height 240mm, Dia.32mm, Up 30 LED & Bottom 92 LED.

SLV-630-J4LDA1 Stand 285*238*40mm, Pole Height 240mm, Dia.32mm, Up 30 LED &

Bottom 92 LED. With
SLV-640-A1 Holder Inner Dia.76mm.



SLV-630-J1/J1A1

SLV-630-J1 Stand 285*238*25mm, Pole Height 240mm, Dia.32mm, No Light Source.

SLV-630-J1A1 Stand 285*238*25mm, Pole Height 240mm, Dia.32mm, No Light Source. With SLV-640-A1 Holder Inner Dia.76mm.

SLV-630 Series Pole Large Stand Models



SLV-630-B3LS/B3LSA1

SLV-630-B3LS Stand 380x260x20mm, Pole Height 240mm, Dia.32mm, Double Pipe LED.

SLV-630-B3LSA1 Stand 380x260x20mm, Pole Height 240mm, Dia.32mm, Double Pipe LED. With SLV-640-A1 Holder Inner Dia.76mm.



SLV-630-B3/B3A1

SLV-630-B3 Stand 380x260x20mm, Pole Height 240mm, Dia.32mm, No Light Source.

SLV-630-B3A1 Stand 380x260x20mm, Pole Height 240mm, Dia.32mm, No Light Source. With A54.3640-A1 Holder Inner Dia.76mm.



SLV-630-B3B/B3BA1/B3BA5

SLV-630-B3B Black Stand 380x260x20mm, Pole Height 240mm, Dia.32mm, No Light Source.

SLV-630-B3BA1 Black Stand 380x260x20mm, Pole Height 240mm, Dia.32mm, No Light Source. With SLV-640-A1 Holder Inner Dia.76mm.

SLV-630-B3BA5 Black Stand 380x260x20mm, Pole Height 240mm, Dia.32mm, No Light Source. With SLV-640-A5 Holder Inner Dia.76mm.

SLV-630 Series Track Stand Models



SLV-630-B4



SLV-630-B7

SLV-630-B4/B7

SLV-630-B4 Track Stand 200x255x22mm, Track Height 300mm, Dia.76mm, Focusing Range 106mm, No Light Source.

SLV-630-B7 Track Stand 205x275x22mm, Track Height 300mm, Dia.76mm, Focusing Range 106mm, No Light Source.



SLV-630-L



SLV-630-2L

SLV-630-L/2L

SLV-630-L Stand 255*200*28mm, Track Height 276mm, Dia.76mm, Focusing Range 78mm, No Light Source.

SLV-630-2L Stand 255*200*40mm, Track Height 276mm, Dia.76mm, Focusing Range 74mm. LED Light.



SLV-630-B9LS

SLV-630-B9L

SLV-630-B9LD

SLV-630-B9LS/B9L/B9LD

SLV-630-B9LS Stand 205x275x40mm, Track Height 300mm, Dia.76mm, Focusing Range 106mm, Double Pipe LED.

SLV-630-B9L Stand 205*275*40mm, Track Height 300mm, Dia.76mm, Focusing Range 106mm, Up/Bottom 3W LED.

SLV-630-B9LD Stand 205*275*40mm, Track Height 300mm, Dia.76mm, Focusing Range 106mm, Up 30 LED & Bottom 92 LED

SLV-630 Series Track Large Stand Models



SLV-630-B5

SLV-630-B5

SLV-630-B5 Stand 380x260x20mm, Track Height 300mm, Dia.76mm, Focusing Range 106mm, No Light Source.

SLV-630 Series Track Fan Stand Models



SLV-630-J2

SLV-630-J2

SLV-630-J2 Stand 285*238*25mm, Track Height 300mm, Dia.76mm, Focusing Range 106mm, No Light Source.

SLV-630 Series Single Arm Stand Models



SLV-630-STL1



SLV-630-STL1A1

SLV-630-STL1/STL1A1

SLV-630-STL1 Stand 230*230*40mm, Pole Height 380mm, Lateral movement 235mm, Dia.32mm, No Light Source.

SLV-630-STL1A1 Stand 230*230*40mm, Pole Height 380mm, Lateral movement 235mm, Dia.32mm, No Light Source.
With SLV-640-A1 Holder Inner Dia.76mm



SLV-630-STL3/STL3A3

SLV-630-STL3 Stand 250*240*12mm, Pole Height 350mm, Lateral movement 400mm, Dia.38mm, No Light Source.

SLV-630-STL3A3 Stand 250*240*12mm, Pole Height 350mm, Lateral movement 400mm, Dia.38mm, No Light Source. With SLV-640-A3 Holder Inner Dia.76mm



SLV-630-STL4/STL4A3

SLV-630-STL4 Stand 260*280*15mm, Pole Height 380mm, Lateral movement 420mm, Dia.37mm, No Light Source.

SLV-630-STL4A3 Stand 260*280*15mm, Pole Height 380mm, Lateral movement 420mm, Dia.37mm, No Light Source. With SLV-640-A3 Holder Inner Dia.76mm

SLV-630 Series Boom Stand Models



SLV-630-STL8/STL8A3

SLV-630-STL8 Stand 380*280*10mm, Pole Height 355mm, Lateral movement 265mm, Dia.32mm, No Light Source.

A54.3630-STL8A3 Stand 380*280*10mm, Pole Height 355mm, Lateral movement 265mm, Dia.32mm, No Light Source.
With A54.3640-A3 Holder Inner Dia.76mm

SLV-630 Series Dual Arm Stand Models



SLV-630-STL2/STL2A1

SLV-630-STL2 Stand 210*260*50mm, Pole Height 380mm, Lateral movement 235mm, Dia.32mm, No Light Source.

SLV-630-STL2A1 Stand 210*260*50mm, Pole Height 380mm, Lateral movement 235mm, Dia.32mm, No Light Source.
With SLV-640-A3 Holder Inner Dia.76mm



SLV-630-STL5

SLV-630-STL5A3

SLV-630-STL5/STL5A3

SLV-630-STL5 Stand 260*280*15mm, Pole Height 380mm, Lateral movement 420mm, Dia.37mm, No Light Source.

SLV-630-STL5A3 Stand 260*280*15mm, Pole Height 380mm, Lateral movement 420mm, Dia.37mm, No Light Source.
With SLV-640-A3 Holder Inner Dia.76mm

SLV-630 Series Swing Arm Stand Models



SLV-630-STL6A

SLV-630-STL6AA3

SLV-630-STL6A/STL6AA3

SLV-630-STL6A Stand 98*80*116mm, vertical Height 400mm, Horizontal distance 770mm, Max clamping 50mm, No Light Source.

SLV-630-STL6AA3 Stand 98*80*116mm, vertical Height 400mm, Horizontal distance 770mm, Max clamping 50mm, No Light Source. With SLV-640-A3 Holder Inner Dia.76mm



SLV-630-STL6B

SLV-630-STL6BA3

SLV-630-STL6B/STL6BA3

SLV-630-STL6B Stand 98*80*116mm, vertical Height 400mm, Horizontal distance 770mm, Max clamping 50mm, Dia.25mm, Pole Height 275mm, No Light Source.

SLV-630-STL6BA3 Stand 98*80*116mm, vertical Height 400mm, Horizontal distance 770mm, Max clamping 50mm, Dia.25mm, Pole Height 275mm, No Light Source. With SLV-640-A3 Holder Inner Dia.76mm



SLV-630-STL6C

SLV-630-STL6CA3

SLV-630-STL6C/STL6CA3

SLV-630-STL6C Stand 260*280*15mm, vertical Height 400mm, Horizontal distance 770mm, Dia.25mm, Pole Height 323mm, No Light Source.

SLV-630-STL6CA3 Stand 260*280*15mm, vertical Height 400mm, Horizontal distance 770mm, Dia.25mm, Pole Height 323mm, No Light Source. With SLV-640-A3 Holder Inner Dia.76mm

SLV-640 Series Holder Arm Models



SLV-640-A1
Pole Dia. 32mm,
Holding Dia. 76mm



SLV-640-A3
Pole Dia. 16mm,
Holding Dia. 76mm



SLV-640-A5
Pole Dia. 25mm,
Holding Dia. 50mm



SLV-640-A5F
Pole Dia. 25mm,
Holding Dia. 50mm
With Coaxial Coarse
& Fine Focusing