



LABMICRO

8553545575 / 8050259428

contact@labmicro.in




www.labmicro.in



DISSECTING MICROSCOPE



MODEL No. SLV-100



Shop No.1, Basement floor, Best Apartments, Dr. Arunachalam Road, 5th Stage BEML Layout,
Rajarajeshwari Nagar, Bengaluru - 560 098



Technical Specification

Description: Dissecting Microscope with Brass Fitting is a great choice for students and hobbyists. It features a high-quality optical system with a 4x objective lens, 10x Lens, The LABMICRO microscope is a top-quality product that offers the best in microscopy technology. It is a reliable and durable product that will give you years of service. The microscope comes with a one-year warranty.

MODEL No.	SLV-100
Viewing Head:	Monocular, 90°
Eyepiece:	WF10x & WF 20X
Specimen stage:	Square stage with slide holder along with Black and white plate
Focusing:	Uniaxial coarse focusing mechanism with limit stop to prevent the contact between Eyepiece and specimen.
Illumination:	Mirror / LED Illumination

OPTIONAL: Pointer Attachment, Micrometre/Grid/Grain Attachment, DIGITAL Attachment,