

SG Series Ice Maker



- Use the imported Speed Reducer from Haitec in Italy and Motor from GGM in Korean.
- Evaporator are made of Copper.
- Compact design , Luxurious appearance.
- Full automatic running and easy to use.
- A power-overloading protecting device, water shortage auto-detecting, system and ice storage full control device.
- High efficiency CFC free compressor, economical and low energy consumption.

Applicable area

Various types of experiments, hospitals, schools, laboratories, medical, chemical applications, supermarkets, hotels and other seafood products stores etc.

| Model | SG-20 | SG-30 | SG-40 | SG-50 | SG-70 | SG-85 | SG-100 | SG-130 |
|-----------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|
| Ice making capacity | 20kg | 30kg | 40kg | 50kg | 70kg | 85kg | 100kg | 130kg |
| Product size | 300X493X54 7mm | 300X493X54 7mm | 380X543X72 2mm | 380X543X72 2mm | 548X611X88 3mm | 548X611X88 3mm | 548X611X88 3mm | 635X611X94 5mm |
| Packing size | 377x571x59 5mm | 377x571x59 5mm | 457x621x77 0mm | 457x621x77 0mm | 625x690x93 2mm | 625x690x93 2mm | 625x690x93 2mm | 712x690x99 4mm |
| Weight | 32kg | 32kg | 40kg | 40kg | 57kg | 57kg | 58kg | 70kg |
| Capacity per day | 20kg/day | 30kg/day | 40kg/day | 50kg/day | 70kg/day | 85kg/day | 100kg/day | 130kg/day |
| Power | 220V/50HZ | 220V/50HZ | 220V/50HZ | 220V/50HZ | 220V/50HZ | 220V/50HZ | 220V/50HZ | 220V/50HZ |
| Cooling | R134a | R134a | R134a | R134a | R134a | R134a | R134a | R134a |
| Ice storage quantity | 10kg | 10kg | 15kg | 15kg | 25kg | 25kg | 25kg | 35kg |
| Doors | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| Production per hour | 0.83kg | 1.25kg | 1.66kg | 2.08kg | 2.92kg | 3.54kg | 4.2kg | 5.4kg |
| Ideal operating temp. | 15°C | 15°C | 15°C | 15°C | 15°C | 15°C | 15°C | 15°C |
| Consumption | 180w | 180w | 280w | 280w | 420w | 420w | 420w | 685w |

