

Inspection Microscope

Zoom Stereo Microscope, 0.7~4.5x

- Continuous Zoom 0.7x-4.5x Binocular & Trinocular Head
- 7x-45x, Up To 3.5x-270x With Optional Eyepieces & Auxiliary Objectives
- High Eyepoint, Dia.30mm WF10x/20mm, WF15x, WF20x, WF25x Optional
- Working Distance 100mm, Up To 30~165mm With Auxiliary Objectives
- 35 Kinds Stand Optional, Details See 3630



SLV-45 Series Zoom Stereo Microscope Specification		B	T
Head	Binocular Head, 45° Inclined, 360° Rotatable, Interpupillary Adjustment 54-75mm, Both Eyepiece Tube ±5° Dioptre Adjustable,	●	
	Trinocular Head, 45° Inclined, 360° Rotatable, Interpupillary Adjustment 54-75mm, Both Eyepiece Tube ±5° Dioptre Adjustable, With Eyepiece Adapter Focus Adjustable, Light Split Fixed P50:E50, Simul Focal Design To Show Image On Both Eyepiece & Photo Port At Same Time		●
Eyepiece	WF10x/20mm, High Eyepoint, Dia.30mm		
Zoom Lens	0.7x~4.5x		
Zoom Ratio	1:6.5		
Magnification	7x-45x, Up To 3.5x-270x With Optional Eyepieces & Auxiliary Objectives		
Working Distance	100mm, Up To 30~165mm With Auxiliary Objectives		
Field Of View	0.9-57.1mm (Detail See Optical Data Sheet)		
Focusing	Pole Stand Coarse Focusing Range 50mm, Tension Adjustable Track Stand Coarse Focusing Range 74mm, Tension Adjustable		
Stand	SLV-630-STL1A1 Stand 230*230*40mm, Pole Height 380mm, Lateral movement 235mm, Dia.32mm, No Light Source. With SLV-640-A1 Holder Inner Dia.76mm		

