

ADVANCE ROTARY MICROTOME

Latest Technology in Manual Rotary Microtome with Separate Feed Range of Slicing & Trimming

Microtome's are widely used in Pathology, Histopathology, Textile and other industries. Sectioning using a semi automatic microtome enables study and analysis in Light microscopy, histology, electron microscopy, Spectroscopy & Botanical Microscopy.

Technical Specifications:

1. Coarse Feed Manual Hand Wheel
2. Manual Forward & Reverse System
3. Range of the thickness of slicing: 0.5 - 60 μm
4. Selection of Slicing Range:
 - 0.5----1 μm increment 0.5 μm
 - 1-----10 μm increment 1 μm
 - 10--- 20 μm increment 2 μm
 - 20----60 μm increment 5 μm
5. Selection of Trimming in Two Steps :
 - Step of 10 Microns
 - Step of 30 Microns
6. Specimen Horizontal Feeding: 28 mm
7. Vertical Stroke Length: 60 mm.
8. Knife Index: High Profile Disposable Blades
9. Retractable specimen clamp, lateral knife holder
10. Weight: 24 kg Approx.
11. Working : Manually Operated
12. Delivery time: 25 Days



MODEL No. SLV-LRM

