



SLV-501C



SLV-501-B



SLV-501-A

| SLV-501 Multi-Viewing Biological Microscope | | A | B | C |
|---|--|---|---|---|
| Head | 2 pcs Sliding Binocular Head, 30°Inclined, Interpupillary Distance 55-75mm | ● | - | - |
| | 2 pcs Monocular Head and 1 pc Sliding Binocular Head, 30°Inclined, Interpupillary Distance 55-75mm | - | ● | - |
| | 5 pcs Sliding Binocular Head, 30°Inclined, Interpupillary Distance 55-75mm | - | - | ● |
| Eyepiece | WF10X | ● | ● | ● |
| Objective | Achromatic 4X | ● | ● | ● |
| | Achromatic 10X | ● | ● | ● |
| | Achromatic 40X(S) | ● | ● | ● |
| | Achromatic 100X(S) Oil | ● | ● | ● |
| Nosepiece | Quadruple Revolving Nosepiece | ● | ● | ● |
| Stage | Double Layer Mechanical Stage 140mm x 140mm, Travel Range 75x50mm | ● | ● | ● |
| Condenser | Abbe N.A.=1.25 with Iris Diaphragm & Filter | ● | ● | ● |
| Focusing | Coaxial Coarse & Fine Focusing Knobs, Fine division 0.002mm, Coarse Stroke 37.7mm Per Rotation, Ffine Stroke 0.2mm Per Rotation, Moving Range 28mm | ● | ● | ● |
| Light Source | Base with Built-in Halogen Lamp 24V/100W , with Continue Intensity Control, 110V/220V | - | - | ● |
| | Base with Built-in Halogen Lamp 6V/20W , with Continue Intensity Control, 110V/220V | ● | ● | - |
| Indicator | Demo. Indicator with Illuminator of 6V/3W | - | - | ● |

