

STDUDENT MICROSCOPES



Monocular Microscope			
Head	Monocular Head,45° inclined		
	Monocular Head,45° inclined,360° Rotatable		
Eyeiece	WF10X		
Objective	Achromatic 4X		
	Achromatic 10X		
	Achromatic 40X, Spring		
	Achromatic 100X, Spring		
Nosepiece	Triple Revolving Nosepiece		
Stage	Double Layer Mechanical Stage with Specimen Clips, Size 90mmx90mm		
Condenser	Condenser NA 0.65		
Focusing	Coarse Focusing		
	Coaxial Coarse & Fine Focusing Adjustment Fine Focusing Scale Value 0.002mm		
Light Source	Up & Bottom Led Illumination.		
Optional Accessories			
Eyeiece	WF16X		
	WF20X		
Objective	Achromatic 20X		
	Achromatic 60X(S)		

Monocular Microscope			
Head	Monocular Head,45° inclined		
	Monocular Head,45° inclined,360° Rotatable		•
Eyepiece	WF10X		•
Objective	Achromatic 4X		•
	Achromatic 10X		•
	Achromatic 40X, Spring		•
	Achromatic 100X, Spring		•
Nosepiece	Triple Revolving Nosepiece		•
Stage	Double Layer Mechanical Stage with Specimen Clips, Size 90mmx90mm		•
Condenser	Condenser NA 0.65		•
Focusing	Coarse Focusing		
	Coaxial Coarse & Fine Focusing Adjustment Fine Focusing Scale Value 0.002mm		•
Light Source	Up & Bottom Led Illumination.		•
Optional Accessories			
Eyepiece	WF16X		
	WF20X		
Objective	Achromatic 20X		
	Achromatic 60X(S)		

