



SLV-135i



SLV-135i + Coaxial Illumination



SLV-135i + FL

Parallel Zoom Fluorescent Stereo Microscope, 0.75x-13.5x, 1:18, APO, Upper & Bottom Light

- Parallel Stereo Microscope With FL, PL, DF, BF/DF, APO Functions
- Zoom Range 0.75x-13.5x, Total Magnification 7.5x-135.0x
- Zoom Ratio 18:1, Covers Macro and Microscopy
- Ergonomic Trinocular Head 20° Inclined with Tilting Heads Options
- Plan-Apo Lens for Excellent Image Quality

	Plan Apo 0.5x	Plan Apo 1x	Plan Apo 1.6x	Plan Apo 2x
Maximum NA	0.075	0.15	0.24	0.3
Working Distance	71mm	60mm	30mm	20mm
Correction Ring	–	–	–	3mm Water Depth
Wavelength	380-700nm			

