

BACTERIOLOGICAL INCUBATOR [SMART]

- Double walled inside aluminium/stainless steel 304 qty.
- Outside of Mild steel sheet furnished in epoxy powder coating.
- Insulation by mineral wool of 75mm.



- Door has double glass fitted for observation of specimen. Glass door inside in stainless steel frame may be provided on request [without glass window outside]
- Heating elements are at bottom and sides.
- Two air ventilators on side.

TEMPERATURE

- Micro processed Electronic digital temperature controller cum indicator from few deg.°C above ambient 80°C with $\pm 0.5^\circ\text{C}$ may be increased up to 100°C on request having two windows display SV & PV, auto turning, timer up 999 Mts [selectable] low & high temp. accuracy alarm, cycle, timer with buzzer.
- Perforated shelves adjustable.
- With air circulating fan.
- Work on 220/230 volts.

SIZE INSIDE CHAMBER [ALUMINIUM]

MODEL	W	D	H	Litres.
ISA-1	300	300	300mm	28
ISA-2	355	355	355mm	45
ISA-3	455	455	455mm	95
ISA-4	455	455	605mm	125
ISA-5	605	605	605mm	224
ISA-6	605	455	910mm	252
ISA-7	605	605	910mm	336

SIZE INSIDE CHAMBER [STAINLESS STEEL]

MODEL	W	D	H	Litres.
ISS-1	300	300	300mm	28
ISS-2	355	355	355mm	45
ISS-3	455	455	455mm	95
ISS-4	455	455	605mm	125
ISS-5	605	605	605mm	224
ISS-6	605	455	910mm	252
ISS-7	605	605	910mm	336