

HOT AIR OVEN UNIVERSAL "Super Deluxe Model"



CONSTRUCTION:-

steel 304 qlty.

epoxy powder coating.

,LHS AND RHS. Ribs.

- Double walled inside aluminum/stainless
- Outside made of Mild steel and painted in
- Insulation of 75mm by mineral wool.
- Heating elements are on three sides ' bottom
- With air ventilators on two sides.

TEMPERATURE:-

Cum Indicator from few °C

Above ambient to 250°C ± 1° C.

Mechanical.

- Microprocessed electronic Digital Controller
- Perforated shelves adjustable.
- Work on 220/230 volts.
- Complete with air circulating fan.
- Optional : Automatic Timer Electronic or

SIZE INSIDE CHAMBER [STAINLESS STEEL]

MODEL	W	D	H	Ltrs.
SLV-HAO1	300	300	300mm	28
SLV-HAO2	355	355	355mm	45
SLV-HAO3	455	455	455mm	95
SLV-HAO4	455	455	605mm	125
SLV-HAO5	605	605	605mm	224
SLV-HAO6	605	455	910mm	252
SLV-HAO7	605	605	910mm	336

Option: Microprocessed PID Controller having 4 ramp and 4 soak profile having set mode, alarm parameter , configuration parameter , profile parameter , set point profile , Pv error indicator, parameter for locking etc.

Alarm : It is possible to provide electronic alarm if temperature deviates more than ±1°C to cost extra.

Automatic timer electronic 999 mts. To extra cost.

