



### **Zoom Stereo Microscope, 0.7~4.5x**

- Continuous Zoom 0.7x-4.5x Binocular & Trinocular Head
- 7x-45x, Up To 3.5x-270x With Optional Eyepieces & Auxiliary Objectives
- High Eyepoint, Dia.30mm WF10x/20mm, WF15x, WF20x, WF25x Optional
- Working Distance 100mm, Up To 30~165mm With Auxiliary Objectives
- 35 Kinds Stand Optional, Details See 3630

## Product Details



Eyepiece WF10X/20mm

0.7x-4.5x, Magnification Up To 3.5x-180x With Optional Eyepieces And Auxiliary Objectives



### 3645, 3645N Series Zoom Stereo Microscope Specification

B

T

3645, 3645N Series Zoom Stereo Microscope Specification		B	T
Head	Binocular Head, 45° Inclined, 360° Rotatable, Interpupillary Adjustment 54-75mm, Both Eyepiece Tube $\pm 5^\circ$ Diopter Adjustable,	●	
	Trinocular Head, 45° Inclined, 360° Rotatable, Interpupillary Adjustment 54-75mm, Both Eyepiece Tube $\pm 5^\circ$ Diopter Adjustable, With Eyepiece Adapter Focus Adjustable, Light Split Fixed P50:E50, <b>Simul Focal Design</b> To Show Image On Both Eyepiece & Photo Port At Same Time		●
Eyepiece	WF10x/20mm, High Eyepoint, Dia.30mm		
Zoom Lens	0.7x~4.5x		
Zoom Ratio	1:6.5		
Magnification	7x-45x, Up To 3.5x-270x With Optional Eyepieces & Auxiliary Objectives		
Working Distance	100mm, Up To 30~165mm With Auxiliary Objectives		
Field Of View	0.9-57.1mm (Detail See Optical Data Sheet)		
Focusing	Pole Stand Coarse Focusing Range 50mm, Tension Adjustable Track Stand Coarse Focusing Range 74mm, Tension Adjustable		
Stand	35 Kinds Stand Optional, Details See <b>3630</b> Choose Model According To Cata. No. Like: <b>3645-B1</b> Means Binocular Head With B1 Stand <b>3645-B1T</b> Means Trinocular Head With B1 Stand		



Shop No.1, Basement floor, Best Apartments, Dr. Arunachalam Road, 5th Stage BEML Layout, Rajarajeshwari Nagar, Bengaluru - 560 098

## 3645 Head Optional



## 35 Kinds Stand Optional

35 Kinds Stand Optional, Details See [3630](#)

Choose Model According To Catalogue No. Like:

**3645-B1** Means Binocular Head With B1 Stand

**3645-B1T** Means Trinocular Head With B1 Stand



## 3645 Series Models



3645-B1T



3645-B2T



3645-B2LT



3645-B3T



3645-B3LST



3645-B2LDT



3645-B4T



3645-B5T



3645-B6T



3645-B7T



3645-B8LT



3645-B8LDT



3645-B9LT



3645-B9LDT



3645-B8LST



3645-LT



3645-2LT



3645-B9LST



3645-J1T



3645-J2T



3645-J3LT



3645-J4LT



3645-J4LDT



3645-J3LDT



Shop No.1, Basement floor, Best Apartments, Dr. Arunachalam Road, 5th Stage BEML Layout, Rajarajeshwari Nagar, Bengaluru - 560 098