



## **Zoom Stereo Microscope, 0.8~5x, Trinocular, Square Stand Up/Bottom LED**

- Zoom Stereo 0.8x~5.0x, Total Magnification Up To 4x~200x With Auxiliary Lens & Eyepiece
- Super Wide Field Plan Eyepiece WF10x/22mm Diopter Adjustable, WF15x, WF20x Optional
- Auxiliary Lens 0.5x, 0.7x, 2.0x Optional To Enlarge Magnification and Working Distance
- Large Track Stand With Up/Bottom LED Light In Square & Fan Shape For Choice
- A56.1522 Extra 60 LEDs Ring Light Source Optional For More Illumination

## Product Details



SLV-50-B3



SLV-50-T4



SLV-50-B1

Track Big Fan Stand,  
No Light



SLV-50-B2

Track Big Fan Stand  
Upper/Bottom LED



SLV-50-B3

Track Big Square  
Stand, No Light



SLV-50-B4

Track Big Square  
Stand, Upper/Bottom  
LED



SLV-50-T1



SLV-50-T2



SLV-50-T3



SLV-50-T4



Standard Eyepiece  
WF10x/22



Optional Eyepiece:  
WF10x/22, WF15x/16, WF20x/12 with Diopter Adjustment



Optional LED Ring Light 60pcs



Auxiliary Objectives: 0.5x/220mm; 0.7x/137mm; 2x/45mm



# TECHNICAL SPECIFICATION

| SLV-50 Series Zoom Stereo Microscope                      |   | B1 | T1 | B2 | T2 | B3 | T3 | B4 | T4 |
|---|---|----|----|----|----|----|----|----|----|
| Head  | Binocular Head, 45° Inclined, 360° Rotatable, Interpupillary Range 52mm~75mm  | o  |    | o  |    | o  |    | o  |    |
|   | Trinocular Head, 45° Inclined, 360° Rotatable, Interpupillary Range 52mm~75mm, Light Spilt Ratio E30:P70, With 0.5x C-Mount |    | o  |    | o  |    | o  |    | o  |
| Magnification   | 8x~50x, Up To 4x~200x With Optional Eyepieces & Auxiliary Objective   |    |    |    |    |    |    |    |    |
| Eyepiece  | WF10x/22mm, Diopter Adjustment -8~+5°   |    |    |    |    |    |    |    |    |
| Zoom Lens   | 0.8~5.0x  |    |    |    |    |    |    |    |    |
| Zoom Ratio  | 1:6.5   |    |    |    |    |    |    |    |    |
| Zoom Knob   | Zoom Knob With <b>Click Stop</b> At 1x,2x,3x,4x   |    |    |    |    |    |    |    |    |
| W.D.  | Working Distance 115mm, Up To 45~211mm With Auxiliary Objective   |    |    |    |    |    |    |    |    |
| Focusing  | Coarse Focusing Knob With Tension Adjustable, Focusing Range 100mm  |    |    |    |    |    |    |    |    |
| Working Stage   | Glass Plate Dia.125mm   |    |    |    |    |    |    |    |    |
| Base & Light Source                                       | Big <b>Fan</b> Base Track Stand, No Light   | o  | o  |    |    |    |    |    |    |
|   | Big <b>Fan</b> Base Track Stand, With Up/Bottom LED   |    |    | o  | o  |    |    |    |    |
|   | Big <b>Square</b> Base Track Stand, No Light  |    |    |    |    | o  | o  |    |    |
|   | Big <b>Square</b> Base Track Stand, With Up/Bottom LED  |    |    |    |    |    |    | o  | o  |
| SLV-50 Series Zoom Stereo Microscope Optional Accessories |   |    |    |    |    |    |    |    |    |
| Eyepiece  | WF10x/22mm, Diopter Adjustment -8°~+5°  |    |    |    |    |    |    |    |    |
|   | WF15x/16mm, Diopter Adjustment -8°~+5°  |    |    |    |    |    |    |    |    |
|   | WF20x/12.5mm, Diopter Adjustment -8°~+5°  |    |    |    |    |    |    |    |    |
| Auxiliary Objective                                       | 0.5x/220mm  |    |    |    |    |    |    |    |    |
|   | 0.7x/137mm  |    |    |    |    |    |    |    |    |
|   | 2.0x/45mm   |    |    |    |    |    |    |    |    |
| Light Source  | LED Ring Light 60 LEDs  |    |    |    |    |    |    |    |    |