

MICROSCOPES



SLV-201

LABORATORY MICROSCOPES |



DIGITAL MICROSCOPES |



FLUORESCENCE MICROSCOPES |



INVERTED MICROSCOPES



POLARIZING MICROSCOPES |



PHASE CONTRAST MICROSCOPES |



UPRIGHT & INVERTED METTALURGICAL MICROSCOPES













DIGITAL STEREO ZOOM MICROSCOPES |















MOTORIZED MICROSCOPES |



PCB INSECTION MICROSCOPES

			
<p>Student Microscope SLV-101</p>	<p>Monocular Microscope SLV-105</p>	<p>Binocular Microscope SLV-20i</p>	<p>Binocular Microscope SLV-20iP</p>
			
<p>Binocular Microscope SLV-21</p>	<p>Trinocular Research Microscope SLV-22</p>	<p>Trinocular Infinity Optical Optics SLV-23</p>	<p>Upright Research Microscope SLV-25</p>
			
<p>Intelligent Revolutionary Products, Operating Comfort ergonomics design. SLV-90</p> <ul style="list-style-type: none"> • NIS60 N-PLN Infinity Plan Objective • SW10x/25mm, High Eyepoint, Diopter Adjustable, Dia.30mm • Coded Nosepiece, Sextuple, Backward, For Auto Brightness Adjust • LABO – 91FN Full Auto Motorized Research Level Scientific Laboratory Microscope • Upgrade To Fluorescent Microscope 		<p>RESEARCH LABORATORY MICROSCOPES SLV-30</p> <ul style="list-style-type: none"> • Infinity Semi- apochromatic Flu Objective • Eyepiece WF10X/25mm (High Eyepoint) • Nosepiece Sextuple, Backward • Working Stage Steel wire stage with Rack-less, Size 280x160mm, moving range 80x52mm. • Upgrade To Fluorescent Microscope 	

REMARKS: -

			
<p>DIGITAL BINOCULAR MICROSCOPE - SLV-20i-3.0M</p>	<p>SLV-21D- 10.5" LCD Touch Screen Digital Microscope, 8.0M, Android Pad, Infinity Plan</p>	<p>SLV-23D-10.5" LCD Touch Screen Digital Microscope, 8.0M, Android Pad</p>	<p>SLV-40- EVOS LCD Cell Image Inverted Microscope, Semi-APO, BF+PH, LCD + Camera</p>
			
<p>SLV-21DA-Digital Biological Microscope, 5.0M</p>	<p>SLV-21D-Digital Biological Microscope, 5.0M</p>	<p>SLV-22D- Biological Digital microscope with Analysis software</p>	<p>SLV-80FL- Digital Fluorescence Microscope with Imagistic</p>
			
<p>SLV-45D- LCD Screen Digital Stereo zoom microscope</p>	<p>SLV-68- Stereo zoom microscope with Mirror Reflection and Illumination moveable</p>	<p>SLV-53T- STEREO ZOOM MICROSCOPE WITH LCD Touch Screen Digital Microscope, 8.0MP</p>	<p>SLV-AMSPIRA Calibration Free Motorized Zoom Video Microscope</p>

<ul style="list-style-type: none"> • Digital Microscope 	<ul style="list-style-type: none"> • Microscope for Inspection 	<ul style="list-style-type: none"> • High-Resolution Microscope
<ul style="list-style-type: none"> • USB Microscope 	<ul style="list-style-type: none"> • Microscope for Inspection 	<ul style="list-style-type: none"> • Microscope for PCB/Electronics



SLV-61I RESEARCH FLUORESCENCE MICROSCOPE, SEMAPO, UPGRADE BF+DF+PL+PH+FL+DIC



SLV-22FL2 LED-FLUORESCENCE MICROSCOPE WITH DUAL FILTER



SLV-22FL4 LED-FLUORESCENCE MICROSCOPE WITH FOUR FILTER



SLV-80FL- DIGITAL FLORESCNCE MICROSCOPE



LED FLOURESCENCE MICROSCOPE

INVERTED MICROSCOPE



SLV-40iFL RESEARCH FLUORESCENCE MICROSCOPE, SEMAPO, UPGRADE BF+PL+PH+FL



SLV-40i BIOLOGICAL INVERTED MICROSCOPE



SLV-40iu- INVERTED METALLURGICAL MICROSCOPE, REFLECT LIGHT, BF/DF

METALLURGICAL MICROSCOPE




SLV-41M INVERTED METALLURGICAL MICROSCOPE, BF, DF, DIC, PL, ECO, SEMI-APO



SLV-40M INVERTED METALLURGICAL MICROSCOPE





UPRIGHT METALLURGICAL MICROSCOPE, REFLECT & TRANSMIT, BF/DF

 <p>SLV-21M</p>	<p>Designed for precision analysis of metal microstructures and surface features. Equipped with reflected light optics for examining opaque specimens. Ideal for metallography, heat treatment studies, and failure analysis. Supports high-resolution imaging and documentation for research and quality control. Widely used in industrial labs, foundries, and material science institutions.</p>
--	--

METALLURGICAL MICROSCOPE REFLECT/TRANSMIT BF/DF

POLARIZING MICROSCOPE

		
<p>TRINOCULAR POLARIZING MICROSCOPE, REFLECT & TRANSMIT LIGHT (SLV-25POL SERIES)</p>	<p>SLV-21POL-T POLARIZING MICROSCOPE, TRANSMIT, TRINOCULAR</p>	<p>SLV-22POL-T&R Polarizing Microscope, Reflect Light, Trinocular</p>

PHASE CONTRAST MICROSCOPE

		
---	--	---

PHASE CONTRAST MICROSCOPE, INFINITY PLAN, TURRET DISC CONDENSER

STEREO ZOOM MICROSCOPE

			
---	---	--	---

<p>ZOOM STEREO MICROSCOPE, 0.7X – 4.5X, TRINOCULAR, SQUARE STAND</p>	<p>STEREO MICROSCOPE TURRET TYPE 2X AND 4X</p>	<p>ZOOM STEREO MICROSCOPE, 0.8~5X, TRINOCULAR, SQUARE STAND UP/BOTTOM LED</p>	<p>ZOOM STEREO MICROSCOPE, 0.65~5.3X, 1:8.1, TRINOCULAR, UP/DOWN 3W LED</p>
---	---	--	--



**INSPECTION MICROSCOPE
0.7X TO 4.5X**

**PARALLEL OPTICAL
STEREO ZOOM
MICROSCOPE – 0.8X
TO 8X**

**MOTORIZED FLUORESCENT STEREO ZOOM MICROSCOPE,
0.63~6.3X, AND 0.8X TO 13.5X**

Industrial & Manufacturing Applications
 PCB Inspection & Soldering
 Quality Control
 Assembly & Rework
 Gemology & Jewellery
 Forensics
 Archaeology & Palaeontology Digital Imaging
 & Documentation
 Photography & Videography
 Software Integration and Dental Labs



**MONO STEREO ZOOM
MICROSCOPE, 0.7X-5.0X**

**NEW 10.1-INCH
HDMI MICROSCOPE
PHOTO
RESOLUTION: 24M,
20M,16M, 12M, 8M,
5M, 3M, 2M**


**2D/3D manual rotate zoom
video microscope
3d 360° observe each side of
object with large depth of field in
real time dynamic image**

**Smart 3D microscope with 16X-
160X Optical system. 8mm
Diagonal FOV**
 • 1920 x 1080 Resolution
 • 3.75µm×3.75µm Pixel
 Size
 • 60fps@2MP Frame Rate
 • HDMI Data interface

REMARKS:-

DIGITAL CAMERAS

			
<p>2MP TO 20MP CMOS CAMERA</p>	<p>HDMI & USB CMOS 5MP CAMERA</p>	<p>6 MP TRUECHROME CAMERA</p>	<p>5MP TO 20MP CMOS CAMERA</p>
			
<p>9" LCD DIGITAL CAMERA, USB2.0 +MOUSE MEASURE, 5.0M</p>	<p>10.5" TOUCH SCREEN ANDROID PAD DIGITAL CAMERA, 8.0M</p>		
			
<p>DIGITAL CAMERA 16MP WITH 11.6" WIN10 TABLET</p>	<p>USB3.0 SONY CMOS COOLED DIGITAL CAMERA, FLUORESCENT IMAGE, GS SHUTTER, 1.7MP~45MP</p>	<p>13.3" WIN11 TABLET TOUCH SCREEN DIGITAL CAMERA, 4K HDMI+USB3.0+WIFI, 25M 4K</p>	
			
<p>Smart 3D microscope camera for metallographic or stereo microscopes. 8mm Diagonal FOV</p>	<p>Maximising Speed & Resolution in sCMOS FL 20 Color Cooled CMOS Camera for Low Magnification Fluorescence Imaging.</p> <ul style="list-style-type: none"> • 15.86 mm Diagonal • 5472 (H) x 3648 (V) • 2.4 μm x 2.4 μm • 14 fps @ 20 MP 		



Overhead Stirrers
High Torque / High Speed




WIGGENS
THE MAGIC MOTION



INCUBATOR / OVEN / VACUUM OVEN



CO2 INCUBATOR / SHEKAR / ROLLER



Homogenizer / Flapping Homogenizer



Shekar Incubator

Shakers
Versatility for a variety of sample preparation and bench applications



ULTRASONIC SHAKER / CLEAN



LAB BURNER / COLONY COUNTER



VISION MEASURING SYSTEM		
STAGE	X/Y Moving scope	200*100
	Z	Available space:220mm, Focusing moving:200mm W D . 95mm
	SIZE	230*130 for Glass plate
	W A Y	XYZ three-axis precision V-type cross guide rails! X/Y axes are driven by optical shafts, and the Z axis is driven by an optical shaft.
	System	Optical Scale Resolution: X/Y/Z 0.0005mm Multifunctional Data Processor Capable of Measuring Points, Lines, Arcs, Circles, Angles, etc
	Measurement Accuracy	≤3μm + L/200 Repeatability: ≤3 μm + L/250
	3D (option)	UK Renishaw Probe Set! Probes come with φ1//φ2
	Vision Targeting System	Camera: Japan Sony Color 1/3" CCD 800TVL Camera, Magnification Ratio 30-190X Lens: WM High-Definition Variable Magnification Lens // Continuous Magnification Ratio 0.7-4.5X
	Light system	Surface Light Source and Transmission Light Source: LED Cold Light Source with Long Service Life and Adjustable Brightness (5-Ring 8-Zone Programmable Light Source) Dual Laser Positioning Indication
	Size (L*W*H)	500*700*1100

MULTI VIEWING MICROSCOPE

DECA HEAD MICROSCOPE



MODEL No. 510 - 2

PENTA HEAD MICROSCOPE



MODEL No. 510 - 1

DECA HEAD MICROSCOPE



MODEL No. 510 - 3

TRI-HEAD MICROSCOPE



			
Leap-Count Cell Counter	Microplate Reader	Elisa Microplate Reader	Triple Gradient Block Thermal Cycler
			
Isothermal Fluorescence PCR	Biosafety Cabinet (Class II B2)	BSC-4FA Series NSF Certified Biosafety Cabinet (Class II A2)	Vertical Laminar Flow Cabinet
			
Horizontal Laminar Flow Cabinet	Ducted PP Fume Hood	Vertical Autoclave (Digital) Capacity: 35L, 50L, 75L, 100L,150L Temp: 105~134°C/0.22Mpa	Vertical Pressure Steam Sterilizer Capacity: 28L Temp: 105~134°C/0.22Mpa
			
Tabletop Autoclave 35L / 50L	-86 Degree Laboratory Deep Freezer 348L	Refrigerator	Scientz-10N Series Freeze Dryer

REMARKS: -

HISTOPATHOLOGY EQUIPMENTS

			
<p>Rotary Microtome</p>	<p>Semi-automatic Microtome + Fast Freezing</p>	<p>Auto Nucleic Acid Purification System</p>	<p>Automated Tissue Processor (Mini-type)</p>
			
<p>Cooling Plate</p>	<p>Automated Tissue Processor (Mini-type)</p>	<p>Fully Automated Tissue Stainer</p>	<p>Slide Dryer</p>
			
<p>Vacuum Automated Tissue Processor</p>	<p>KD-HI Slide Dryer</p>	<p>Benchtop pH/Ion/Conductivity/DO Meter</p>	<p>Rotary Evaporator</p>
	<p>Chair & Stool</p> 		
<p>Customized LAB TABLE SETUP / Furniture's /CHAIRS</p>		<p>SS 304 or 306 CABINETS</p>	<p>Ice Maker</p>

SPECIALIZED METALLOGRAPHY INSTRUMENTS



Metallurgical Cutting Machine Manufacturer



Mounting Press



Full Rockwell Hardness Tester



Alpha-660/661 Dual-Disc Automatic Grinder and Polisher




12MHBS-3000ZD Intelligent automatic Brinell hardness tester



HV-50A Automatic turret Vickers hardness tester



Cut Off Wheels



Metal Bond Wafering Blades




Cutting Coolant



Cutting

Accessories

Suppliers Hot Mounting Consumables







Research 3D Full Auto Super
Depth of Field Upright
Metallurgical Microscope,
T/RBF+DF+PL+DIC, Semi-APO



All-in-One Industry Microscope, Motorized XYZ + Motorized Nosepiece + Marco Camera



Motorized Fluorescent Microscope,
BF/FL, XYZ Motorized



			
<p>LABORATORY HOMOGENIZER</p>	<p>MAGNETIC STIRRER COMES WITH A HEATING MODULE</p>	<p>HEAT-COLLECTING MAGNETIC STIRRER</p>	<p>CONTINUOUS PIPETTE TIP</p>
			
<p>WATER QUALITY ANALYZER</p>	<p>HEATING MAGNETIC STIRRER</p>	<p>INFRARED INOCULATING LOOP STERILIZER</p>	<p>MAGNETIC STIRRER</p>
			
<p>OIL-FREE DIAPHRAGM VACUUM PUMP</p>	<p>THERMOSTATIC SHAKING WATER BATH</p>	<p>FOUR POSITIONS MAGNETIC STIRRER</p>	<p>AIRFLOW DRYER</p>
			
<p>DIGITAL BOTTLE TOP LIQUID DISPENSER</p>	<p>WATER BATH</p>	<p>MULTIPLE MAGNETIC STIRRER WITH HEATING</p>	<p>DRY BATH</p>


REMARKS:-



ANALYTICAL BALANCE



ELECTRONIC BALANCE



DUAL RANGE ANALYTICAL BALANCE



Overhead Stirrer



CGF CENTRIFUGE



MICROPLATE CENTRIFUGE



LOW/HIGH SPEED CENTRIFUGE



MINI CENTRIFUGE - MC-5PRO/MG7PRO



MINI VACUUM PUMP



LARGE CAPACITY MAGNETIC STIRRER 10L/15L/30L/50L



MAGNETIC STIRRER



HOTPLATE Multi Plate Vortex Mixer



MAGNETIC STIRRER



ROLLER MIXER



ROTATING MIXER SHAKER



Multi Tube Vortex Mixer



LINEAR/ CIRCULAR SHAKER




DRY BATH SHAKING INCUBATOR

REMARKS:-

Professional & Friendly Message

"We offer a wide range of high-quality microscopes and laboratory instruments, along with reliable after-sales service. Whether you're making a new purchase or need support with existing equipment, our team is here to assist you every step of the way. Your satisfaction and trust are our top priorities."



LAB EQUIPMENT REPAIR & SERVICES

Providing Quick and Reliable Service
Beyond Manufacturer Support



MICROSCOPE SERVICING
Routine Servicing - Cleaning, Maintenance & Repairs

Formal Business Message

"We specialize in the supply and servicing of premium microscopes and laboratory instruments. Our commitment to quality products and professional service ensures long-term satisfaction and trust. Contact us today for sales inquiries or service support."



SLV TECHNO LABMICRO
INDIA PVT LTD

<p>SLV TECHNO No.1, Basement Floor, Best Apartment, 5th Stage BEML Layout, Raj Rajeshwari Nagar, Bangalore-560098</p>	<p>SLV TECHNO Plots no 271 ,A.N.S Nagar, Annexe 1, Palavedu village Pakkam - 602024. Tiruvallur district</p>
<p>SLV TECHNO T.C.20/2716 (1), Nedumcaud Road, Karamana, Thiruvananthapuram 695002 Kerala</p>	<p>SLV TECHNO #101, Eunice Building, Purna Bazar East Dimapur 797116 Nagaland</p>