

STEREO ZOOM MICROSCOPE



Stereo Zoom Microscope with 0.65-6.5x Magnification and Ratio 1:10

- Continuous Zoom 0.65x-6.5x Binocular & Trinocular Head, Zoom Ratio 1:10
- Standard 6.5x-65x, Up To 1.95x-325x With Optional Eyepieces & Auxiliary Objectives
- High Eyepoint, Dia.30mm WF10x/22mm, Optional WF15x, WF20x, WF25x
- Working Distance 110mm, Up To 30~287mm With Auxiliary Objectives
- 35 Kinds Stand Optional,



Stand- B6



SLV-65B8LT



Stand-B8LS




SLV-65B3T



TECHNICAL SPECIFICATION

Optical System	Greenough Optical System
Head	Binocular / Trinocular Head, Interpupillary Distance 55~75mm
Magnification	0.65x~6.5x Continuous Zoom Body, Total Magnification 1.95x~325x (With Optional Eyepiece & Lens)
Zoom Ratio	1:10
Eyepiece	WF10x/22mm, High Eyepoint, Diopter Adjustable, Dia.30mm
Max Field of View	Φ36 mm
Working Distance	110mm, Optional 30~287mm
Bracket Holder Hole	Inner Dia. Φ76mm
Model	Trinocular Head, Light Split Fixed 4/6 B6 Pole Stand, Size 210*274*38mm, Pole Dia.32mm, Height 240mm, No Light Source, 0.5x C-Mount
	Trinocular Head, Light Split Fixed 4/6, B8L Pole Stand, Size 210*274*38mm, Pole Dia.32mm, Height 240mm, Up/Bottom 3W LED, 0.5x C-Mount
	Trinocular Head, Light Split Fixed 4/6, B8LS Pole Stand, Size 210*274*38mm, Pole Dia.32mm, Height 240mm, Up Dual Pipe LED Brightness Separately Adjustable, Bottom 3W LED. 0.5x C-Mount
	Trinocular Head, Light Split Fixed 4/6, B3 Pole Stand, Size 380*260*20mm, Pole Dia.32mm, Height 240mm, No Light Source. 0.5x C-Mount



SLV-65T Microscope with B3T Stand along with Digital WIFI Eyepiece Camera



SLV-65 -B3T



SLV-65T STEREO ZOOM MICROSCOPE OPTICAL DATA SHEET

EYEPIECE	Auxiliary Objective Lens											
	0.3x		0.5x		0.7x		1.0x (Standard)		1.5x		2.0x	
	Working Distance		Working Distance		Working Distance		Working Distance		Working Distance		Working Distance	
	287mm		165mm		120mm		110mm		45mm		30mm	
	Mag.	View(mm)	Mag.	View(mm)	Mag.	View(mm)	Mag.	View(mm)	Mag.	View(mm)	Mag.	View(mm)
Wf10x/22mm	1.95x	107.3	3.25x	64.4	4.55x	46	6.5x	32.2	9.75x	21.5	13x	16.1
	19.5x	11.3	32.5x	6.8	45.5x	11.2	65x	3.4	97.5x	2.3	130x	1.7
WF15x/16mm	2.92x	78.03	4.87x	46.83	6.82x	33.45	9.75x	23.42	14.62x	15.64	19.5x	11.71
	29.2x	8.21	48.75x	4.95	68.25x	8.15	97.5x	2.47	146.25x	1.67	195x	1.24
WF20x/12mm	3.9x	58.52	6.5x	35.13	9.1x	25.09	13x	17.56	19.5x	11.73	26x	8.78
	39x	6.16	65x	3.71	91x	6.11	130x	1.85	195x	1.25	260x	0.93
WF25x/9mm	4.87x	43.89	8.12x	26.34	11.37x	18.82	16.25x	13.17	24.37x	8.79	32.5x	6.59
	48.75x	4.62	81.25x	2.78	113.75x	4.58	162.5x	1.39	243.75x	0.94	325x	0.69

DIFFERENT TYPES OF MICROSCOPES STANDS AVAILABLE

w



Shop No.1, Basement floor, Best Apartments, Dr. Arunachalam Road, 5th Stage BEML Layout, Rajarajeshwari Nagar, Bengaluru - 560 098