



LABMICRO

8553545575 / 8050259428 

contact@labmicro.in 

www.labmicro.in 

STEREO MICROSCOPE



**MODEL
SLV-12**



**Shop No.1, Basement floor, Best Apartments, Dr. Arunachalam Road, 5th Stage BEML Layout,
Rajarajeshwari Nagar, Bengaluru - 560 098**

Product Details



01 BINOCULAR HEAD

45° inclined, with WF10x eyepiece, left & right eye tube diopter adjustable +/-5°.

02 OBJECTIVE

Two different models for choosing (1x, 2x, 3x & 1x, 2x, 4x) according to your needs.



03 TRACK STAND


With portable handle.

04 LIGHT SOURCE

Top 1W LED.

Bottom 7x20mA LED, Brightness Adjustable.

Working Stage: Black/white plastic plate dia. 57mm 1pc, frosted glass plate 1 pc.





05

POWER SUPPLY

AC 100~240V, with 3xAA rechargeable batteries built-in base, for order shipped by air, the batteries will not be included.

Optional Accessories

Optional Accessories	
Eyepiece	WF5x/22mm
	WF10x/20mm
	WF15x/15mm
	WF20x/10mm

WF5X/22mm

WF10X/20mm



WF15X/15mm

WF20X/10mm



Product Display



SLV-12-TA



SLV-12-TB

Product Specification

Items	SLV-12-TA	SLV-12-TB
Head	Binocular Head, 45° Inclined, Interpupillary Distance 48~75mm	
Eyepiece	WF10x/20mm	
Objective	1x,2x,3x or 1x,2x,4x (Optional)	
Focusing	Coarse Focusing, Range 40mm	
Working Distance	57mm	
Stand	Track Stand + Touch Switch	Track Stand + Dial Switch
	With Dia 70mm Glass Stage and White/Black Stage	
Light Source	Top 1W LED, Bottom 7x20mA LED Light Brightness Adjustable	
Power	AC 100~240V Wide Range Power, With 3xAA Rechargeable Batteries Built-in Base. (For Order Shipped by Air, The Batteries Will Not Be Included)	